


ArthroFree[®]

The Wireless Camera for Arthroscopy and Endoscopy



 **lazurite[®]** | agile technologies
for optimal care



Wire-Free Simplicity, for Efficiency and Safety¹

Surgeons asked. Nurses dreamed. We delivered—with projected cost-per-case savings of 20%.²


Efficiency at every angle



Time

The ArthroFree System:

- Offers drop-in compatibility with contemporary OR technologies.
- De-clutters the sterile field.
- Reduces set-up and take-down times.³



Energy

One is the magic number

- One sterilization tray for ArthroFree components
- One sterilization process, requiring one chemical treatment

Low wattage, smart consumption

- Conventional light box (90W to 400W): in use all day long
- ArthroFree light source (5W), with rechargeable battery: in use from white balance to end of procedure



Value

The efficiency of simplicity⁴

- Fewer components on the surgical tower
- Fewer items in the sterile field
- No cables to prep, untangle, or replace

Safety, cool and simple



The ArthroFree System eliminates camera-cabled-related:⁵

- Patient burns.
- OR fires.
- Staff trips and falls.

The ArthroFree Advantage

Wireless technology for agility and security

- Surgeon-approved ergonomics and more freedom of movement⁶
- Proprietary wireless connection for patient-data security
 - HIPAA-compliant, point-to-point encryption
 - Not dependent on Bluetooth or Wi-Fi

Low-heat efficiency for patient safety⁷


- Peak recorded temperature of fiber-optic wire 194.5 °C
- Drape ignition flashpoint 121.1 °C
- ISO standard 60601-1: Allowable max. temperature for device contact with patient 48.0 °C
- **ArthroFree's Meridiem® light engine, maximum temperature 43.0 °C**

The five-step surgery prep


- 01** Remove the sterile components from their tray
- 02** Attach an endoscope and slide the Meridiem light engine onto the scope's ACMI input
- 03** Insert a charged battery to power on the system
- 04** Adjust the white balance
- 05** Experience the comfort and freedom of ArthroFree

Less is more ...
with ArthroFree

Conventional System



The ArthroFree System



Technical Specifications

The ArthroFree OR Kit includes: Camera Heads (3), Receiver (1), Batteries (15), Battery Charger (1), and Sterilization Trays (3). All Class II components, except for the Battery Charger. (Images not to scale.)



Camera Head

Imaging	1/2.7" CMOS HD
Sterilization	STERRAD®
Dimensions (handle)	13.8 x 5.0 x 4.2 cm (6.4 cm at antenna)
Weight, with battery	438 g
Product #	AF0001-101-11



Receiver

Imaging	1080p @ 60 fps HD video
Trigger ports	2x 3.5 mm mono (1/8") female
Video outputs	2x DVI
Dimensions (front)	16.4 x 11.3 x 4.2 cm
Weight	294 g
Product #	AF0002-101-10



Battery

Output electrical ratings	3.6 V 2930 mAh, 10.55 Wh
Sterilization	STERRAD®
Dimensions	11.2 x 3.8 x 2.4 cm
Weight	83 g
Product #	AF0003-101-10



Battery Charger

Output electrical ratings	12 vdc, 1 A, 5 bay
Dimensions, with batteries	35.7 x 14.9 x 10.6 cm
Weight, without batteries	800 g
Product #	AF0004-101-10



Sterilization Tray

Dimensions	32.6 x 18.4 x 6.8 cm
Weight	1400 g
Product #	AF0005-101-10

Notes

¹ 82% of respondents (n=84) indicate that they expect the ArthroFree System to improve safety. 84% of respondents (n=85) indicate that they expect the ArthroFree System to improve efficiency. Data on file. Lazurite LLC. Human Factors Summative Usability Studies. 2021-2022.

² This analysis includes amortized capital, sterilization costs, disposable/replaceable items, and service-contract costs. The projected 20% is based upon 30 arthroscopic procedures per week. Lazurite LLC. 22-015R. Cost-Per-Case Analysis of Surgical Camera Systems: The Wireless ArthroFree™ System vs. Conventional Wired Systems vs. Disposable Surgical Camera Systems. 2022.

³ The setup, takedown, and processing times of the ArthroFree System were significantly lower than respective times of a conventional camera system. Lazurite LLC. 21-073R. Simulated OR Wireless versus Wired Arthroscopic Camera Set-Up and Take Down: Summary Report. 2022.

⁴ The ArthroFree System typically replaces two box components on the surgical tower, light wires, and cable guides with a receiver and a battery.

⁵ Bullet 1: Maximum heat production of the ArthroFree wireless camera: 43°C. Select legacy cameras exceed 190°C. Data on file. Lazurite LLC. 21-046R. External Temperature Validation. 2021. Bullet 2: Cowles CE Jr, Culp WC Jr. Prevention of and response to surgical fires. *BJA Education*. 2019; 19(8):261-266.

⁶ Data on file. Lazurite LLC. 21-043R. Human Factors Summative Usability Study of the ArthroFree Wireless Surgical Camera System. 2021.

⁷ Bullet 1: Spradling K, Uribe B, Okhunov Z, et al. Evaluation of ignition and burn risk associated with contemporary fiberoptic and distal sensor endoscopic technology. *Journal of Endourology*. 2015; 29(9): 1076-82. Bullet 2: 3M. "Article Information Sheet. 3M™ loban™ 2 Antimicrobial Incise Drapes 6035, 6040, 6048, 6050, 6051, 6640EU, 6648EU, 6650EU, 6651EU." Document Group: 08-3227-9. Version: 1.00. October 13, 2021. <https://tinyurl.com/3Msafetydata>. Bullet 3: ISO standard at 48°C for contact of 1m 1s to 9m 59s. SL Power Electronics Corp. "Considerations for Maximum Allowable Temperature as per Safety Standards IEC 60601-1, IEC 60950-1 and IEC 61010-1." AN-G012. Application Note V0719. 2019. Bullet 4: Data on file. Lazurite LLC. 21-046R. External Temperature Validation. 2021.

Learn more about the advantages of ArthroFree.

Contact us at [1.833.214.2324](tel:1.833.214.2324) or sales@lazurite.co.

lazurite.co/arthrofree

22006D/20231.0

© 2023 Lazurite LLC. All Rights Reserved.

Lazurite, ArthroFree, Meridiem, and the Lazurite logo are registered trademarks of Lazurite.

All other brand and product names are trademarks or registered trademarks of their respective holders.