

Cut the cord

Elevate your endoscopic efficiency

ArthroFree®

The Cordless Endoscopic
Camera for ENT

With **ArthroFree**, cords and bulky light and power boxes are a thing of the past. Experience a new world of flexibility, not only with more freedom of movement, but the ability to expand your endoscopic setup.

Uninterrupted **focus** on what truly matters - patients

- ✓ **Simple** - Eliminate camera-cable management & interruptions
- ✓ **Clear** - Delivers a crystal clear, crisp surgical image
- ✓ **Focus** - Fix your eyes on the screen—and keep them there

Refined workflow, refined **efficiency**

- ✓ **Quick setup** - Attach your scope, insert the battery, and begin your procedure
- ✓ **Increased patient comfort** - No draped cables over your patients allows for an unobstructed, more relaxed, experience
- ✓ **Compatibility** - Designed to integrate seamlessly with both rigid and flexible scopes



SCHEDULE A DEMO TODAY

Chuck Veselits
Director of Business Development
Lazurite Inc.
(216) 870-9551
c.veselits@lazurite.co

Rachel Gajda, COPM-C
National Market Manager
Tungsten Medical (Nat'l Distributor)
(734) 718-9297
rachel@tungstenmed.com

ArthroFree[®]

Endoscopy Tower Kit MTV0062

- 1 ArthroFree Wireless Camera Head**
1080p HD camera with a Meridiem[®] Light Engine. Paired with optimized image processing firmware and white balance algorithm, it ensures a crisp, clear surgical image.
- 2 Wireless Receiver**
Uses encrypted ultra-wide band (UWB) wireless connection, avoiding reliance on customer IT systems. With 100ms video latency at 30fps, you can expect smooth and precise real-time visualization.
- 3 Charger and Batteries**
Each battery delivers 60 min of power. The 5-bay charger ensures your camera is always ready for use, keeping your busy clinic days uninterrupted.
- 4 8MP Clinical Review Display Monitor**
27-inch 4K 8MP IPS display with 178° of viewing.
- 5 Tryten Tower**
Mobile cart option that allows the camera system to travel between exam rooms, optimizing your clinical workflow.



MTR133 Recorder

Next generation DICOM, AI and 4K ready recorder with automated workflow. It's 13.3" touchscreen is an ideal solution for high-quality video capture for comprehensive documentation.

